# ENGIpore Spine

CE 0373

---

Italian IT

English EN

French FR

German DE

Spanish ES

Portuguese PT

Turkish TR

Dutch NL

Bulgarian BG

Slovak SK

Czech CZ

Slovenian SLO

---

|   |  |
|---|--|
|    | <p>Fabbricante<br/> Manufacturer<br/> Fabricant<br/> Hersteller<br/> Fabricante<br/> Fabricante<br/> Üretici<br/> Fabrikant<br/> Производител<br/> Výrobca<br/> Výrobce<br/> Proizvajalec</p>  |
|    | <p>Non risterilizzare<br/> Do not resterilize<br/> Ne pas restériliser<br/> Nicht reesterilisieren<br/> No volver a esterilizar<br/> Não submeter à reesterilização<br/> Sterilize etmeyiniz<br/> Niet opnieuw steriliseren<br/> Да не се стерилизира повторно<br/> Nesterilizujte opakovane<br/> Ne re-sterilizovat<br/> Ne sterilizirajte znova</p>  |
|    | <p>Non riutilizzare<br/> Do not re-use<br/> Ne pas réutiliser<br/> Nicht wiederverwenden<br/> No volver a usar<br/> Não reutilizar<br/> Yeniden kullanmayınız<br/> Niet hergebruiken<br/> Да не се използва повторно<br/> Neroužívejte opakovane<br/> Pouze na jednorázové užití<br/> Samo za enkratno uporabo</p>   |
|  | <p>Non utilizzare se la confezione è danneggiata<br/> Do not use if the package is damaged<br/> Ne pas utiliser si l'emballage est endommagé<br/> Nicht benutzen, wenn die Verpackung beschädigt ist<br/> No se debe utilizar si el envase está dañado<br/> Não utilizar se a embalagem estiver danificada<br/> Ambalaj hasarlı ise kullanmayınız.<br/> Niet gebruiken als de verpakking beschadigd is<br/> Да не се употребява, ако опаковката е повредена!<br/> Neroužívejte, ak je obal poškozený<br/> Při poškozeném obalu nelze použít<br/> Izdelka ne uporabljajte, če je embalaža poškodovana</p> |

|   |   |
|---|---|
|    | <p>Sterile, incluso il metodo di sterilizzazione (radiazioni ionizzanti)<br/> Sterile, including the method of sterilization (ionizing radiations)<br/> Stérile, incluant la méthode de stérilisation (radiations ionisantes)<br/> Steril, samt Sterilisationsmethode (ionisierende Strahlungen)<br/> Estéril, incluido el método de esterilización (radiaciones ionizantes)<br/> Estéril, incluindo o método de esterilização (radiações ionizantes)<br/> Steril, sterilizasyon yöntemi (iyonlaştırıcı radyasyon) dahil<br/> Steriel, inclusief sterilisatiemethode (ioniserende straling)<br/> Стерилен, включително метод на стерилизация (ионизиращо лъчение)<br/> Sterilné, metóda sterilizácie (ionizujúce žiarenie)<br/> Sterilní, včetně sterilizační metody (ionizující záření)<br/> Sterilno, vključno z načinom sterilizacije (ionizirajoče sevanje)</p> |
|    | <p>Utilizzare entro<br/> Use by<br/> Date limite d'utilisation<br/> Mindestens haltbar bis<br/> Fecha de caducidad<br/> Usar antes de<br/> Son kullanma tarihi<br/> Te gebruiken tot<br/> Используй до<br/> Dátum expirácie<br/> Datum použitelnosti<br/> Rok uporabnosti</p>   |
|    | <p>Marcatura CE<br/> CE marking<br/> Marquage CE<br/> CE-Kennzeichnung<br/> Marcado CE<br/> Marcação CE<br/> CE işareti<br/> CE-markering<br/> CE маркировка<br/> CE značka notifikovaného organu<br/> CE značka<br/> Oznaka CE</p>   |
|  | <p>Numero di lotto<br/> Batch code<br/> Numéro de lot<br/> Chargennummer<br/> Número de lote<br/> Número de lote<br/> Parti numarası<br/> Lotnummer<br/> Партиден №<br/> Číslo šarže<br/> Číslo šarže<br/> Številka serije</p>  |

|   |  |
|---|--|
|    | <p>Numero di catalogo<br/>         Catalogue number<br/>         Numéro de catalogue<br/>         Katalog-Nummer<br/>         Número catálogo<br/>         Número de catálogo<br/>         Katalog numarası<br/>         Catalogusnummer<br/>         Каталоген №<br/>         Katalógové číslo<br/>         Katalógové číslo<br/>         Kataloška številka</p>  |
|    | <p>Numero di pezzi<br/>         Number of pieces<br/>         Nombre de pièces<br/>         Stückzahl<br/>         Número de piezas<br/>         Número de unidades<br/>         Adet sayısı<br/>         Aantal stuks<br/>         Бройки от продукта<br/>         Počet kusov v balení<br/>         Počet kusů<br/>         Število kosov</p>  |
|    | <p>Attenzione, vedere le istruzioni per l'uso<br/>         Warning: see instructions for use<br/>         Attention, voir les instructions pour l'emploi<br/>         Achtung, siehe die Gebrauchshinweise<br/>         Cuidado, consultar las instrucciones de uso<br/>         Atenção, consultar as instruções de uso<br/>         Dikkat, kullanım talimatlarına bakınız<br/>         Let op, zie de gebruiksaanwijzing<br/>         Внимание, прочетете инструкцията за употреба!<br/>         Pozor, přeštudujte si priložený návod na použití<br/>         Pozor, přečtěte si návod k použití<br/>         Opozorilo: upoštevajte navodila za uporabo</p> |
|  | <p>Consultare le istruzioni operative<br/>         Consult the operating instructions<br/>         Voir les instructions pour l'utilisation<br/>         Bedienungsanleitung beachten<br/>         Consultar las instrucciones operativas<br/>         Consultar as instruções operacionais<br/>         İşletim talimatlarına başvurunuz<br/>         Lees de gebruiksaanwijzing<br/>         Консултирайте се с инструкцията за употреба!<br/>         Dodržujte návod na použití<br/>         Prostudujte si pokyny<br/>         Preverite navodila za uporabo</p>  |

# ENGPore Spine

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# ENGPore Spine

## High porosity bioceramic bone substitutes

### DESCRIPTION

Hydroxyapatite is found in human bones in percentages ranging between 65% and 80% and it represents the most important constituent of the bone structure. ENGPore bone substitutes are made of biomimetic porous hydroxyapatite, which is very similar in its chemical composition and microstructure to the mineral component of vertebrate bones and teeth. These medical devices are characterized by a very high porosity, which can reach 90% relative to the total volume.

Interconnected porosity is mainly composed of macro- and micro-pores. The size of macro-pores ranges mainly between 200 and 500  $\mu\text{m}$ , while the size of interconnection micro-pores ranges between 80-200  $\mu\text{m}$ . Moreover, these devices feature an intergranular microporosity (<10 micron), which facilitates the absorption of physiological fluids. The morphological features of ENGPore bone substitutes promote cell migration and adhesion for mineral matrix synthesis, thereby promoting rapid bone regeneration.

### INDICATIONS for USE

ENGPore bone substitutes can be used to recover bone substance loss, in particular:

- in presence of bone defects due to trauma, degenerative diseases of the skeletal system or surgical excisions;
- spinal arthrodesis procedures (interbody fusion and posterolateral fusion), and for spondylolisthesis;
- in all pathological conditions in spinal surgery that require the generation of new bone tissue.

ENGPore can also be used:

- in combination with other Finceramica bone substitutes;
- in combination with autologous bone;
- mixed with marrow blood, platelet gel and growth factors;
- mixed with haemostatic agents and biological glues;
- mixed with antibiotics.

### AVAILABLE SIZES

ENGPore is available in pre-shaped forms and in chips.

## INSTRUCTIONS for USE

The following instructions are to be considered as recommendations for the surgical use of ENGIpore; they are not intended for modifying or replacing standard procedures for the treatment of bone defects. Surgical procedures involving the use of bone substitutes may have variable outcomes. The factors to be taken into account when choosing a bone substitute and the surgical technique are as follows: age and bone quality of the patient, location of the defect, expected loading condition, level and quality of tissue vascularity near the defect, possibility for the substitute to be in direct contact with the vital bone tissue, chances to ensure the complete filling of the bone defect, possibility to obtain correct and sufficient primary stabilisation of the bone substitute, chances to obtain complete coverage of the bone substitute.

### **Pre-operative treatment**

It is advisable to carry out pre-operative radiographs to evaluate the conditions of the bone segment and the size of the defect to be treated as well as to choose the type and amount of product to apply and the fixation methods and materials. In the event the surgeon decides to proceed with surgery, the patient must undergo antibiotic prophylaxis. An alternative prophylaxis is essential to be used in patients allergic to penicillins.

It is important to carefully verify the absence of infections or inflammations prior to surgery.

### **Intra-operative aspects**

Use of ENGIpore bone substitutes only for surgeries performed in sterile environments, adopting every precaution. The product must be opened only in sterile environments.

ENGIpore does not interfere with normal surgical procedures. After careful cleaning of the bone surface, the product can be directly applied onto the bone area to be treated.

As for the pre-shaped formats, ensure maximum contact surface at the interface between ENGIpore and the recipient bone in order to guarantee efficient osteointegration.

Do not apply too much pressure when placing the implant; due to its porosity, incorrect handling may cause microfractures in the product.

Any dispersion of material outside the implant site should be carefully avoided, but the amount of product used must be enough to fill the entire bone defect. Excess material must be carefully removed making sure not to dislocate the bone substitute placed in the defect.

Once the filling is completed, close the surgical site following the standard procedure.

ENGIpore does not have sufficient mechanical resistance to bear direct load in the immediate post op phase. If applied in areas that require immediate stabilisation and /or loading, it is important to use fixation means or synthetic systems.

### **Post-operative treatment**

The patient should undergo a suitable post-operative antibiotic therapy according to the surgical procedures. The surgeon must give the patient all the instructions related to the specific application to guarantee a proper post-operative course.

It is recommended to verify the newly formed bone tissue radiographically before subjecting the area to direct loading.

## **CONTRAINDICATIONS**

ENGIpore should not be used in presence of inflammations and infections in the surgical site to be treated.

As for neurosurgical and spine applications, ENGIpore must not be used in presence of dura mater injuries.

Moreover, it should not be used in patients hypersensitive to calcium phosphate salts or affected by: severe diabetes and demineralisation, serious impairment or liver and renal diseases, hypercalcaemia, thyroid and parathyroid complications, altered coagulation, infections due to previous surgeries, metabolic disorders of any kind, nor in patients undergoing steroid and anticoagulant therapies.

The product must not be implanted in areas where complete coverage with soft tissue is not possible.

No side effects, due to the use of ENGIpore in presence of cancer have been reported. However, patients undergoing bone filler treatment following tumour removal must be monitored.

## **WARNINGS**

ENGIpore does not have sufficient mechanical resistance to bear direct loads before regeneration of the newly formed tissue. If applied in areas that require immediate stabilisation and/or loading, it is important to use fixation means or synthetic systems.

ENGIpore is intended to be used only by qualified and specialised medical personnel.

Do not use the product after its expiry date.

Do not use the product if the packaging is damaged.

ENGLpore is a single-use product. Do not use any remains for other operations or other patients but dispose of it as contaminated hospital waste according to the regulations in force at the hospital facility. The product must not be re-sterilised.

The manufacturer is not liable for risks deriving from processes to recondition the device; reconditioning of the device implies the loss of its features and functions.

The package contains an additional label displaying product tracking information.

the implant must be tracked with the label which must be applied to the patient's medical records.

### **SIDE EFFECTS**

No side effects have been reported following the use of this product according to the correct indications.

### **STERILIZATION**

All ENGLpore formats are supplied sterile (sterilized by 25 kGy gamma rays). These products are for single-use and must not be re-sterilised.

Its reuse, in whole or in part, may involve the risk of cross-contamination and the risk of infection at the implant site.

### **STORAGE**

Keep the product in a cool and dry place, protected from light and from heat sources (+10 °C / +40 °C).

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**FIN-CERAMICA FAENZA S.p.A.**