MaioRegen

€ 0373

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Chondral and Osteochondral Substitutes

DESCRIPTION

MaioRegen is a line of reabsorbable implantable medical devices, intended for treating lesions of the cartilage which may affect the subchondral bone structure.

MaioRegen medical devices have a multi-layer structure composed of a layer of collagen fibre (type I, of equine origin) and one or two composite layers of collagen fibre (type I, of equine origin) and magnesium-enriched hydroxyapatite crystals; these are biomimetic materials which are physiologically present in the organism, and which therefore do not alter the biochemical equilibrium of the joint.

All MaioRegen devices have two distinct surfaces: the part composed of collagen is smooth, while the composite part is rough.

MaioRegen devices take the form of a porous composite matrix which mimics the cartilage or the entire bone/cartilage system.

The high porosity and hydrophilicity of MaioRegen devices promotes rapid absorption of biological fluids and the associated cellular components with which they come into contact, thus promoting their adaptation to the implant site.

Thanks to their special configuration and biomimetic qualities, these devices promote the regeneration of cartilaginous and osteo-cartilaginous tissue, and are reabsorbed as they remodel the newly formed tissue.

INDICATIONS for USE

MaioRegen Prime is indicated for the treatment of single or multiple grade IV (Outerbridge classification) osteochondral defects with severely compromised bone.

Such lesions may be traumatic, post-traumatic, degenerative or due to osteochondritis dissecans (OCD).

MaioRegen Prime is also indicated for early OA lesions, as defined (Luyten et al., Knee Surg Sports Traumatol Arthrosc., 2012), of grade I and II according to their classification by Kellgren and Lawrence, in the absence of osteophytes.

Maioregen Slim is indicated for the treatment of single or multiple grade III or IV (Outerbridge classification) deep chondral and osteochondral lesions with slightly compromised bone.

The lesions may be traumatic, post-traumatic or degenerative.

Maioregen Slim is also indicated for early OA lesions, as defined (Luyten et al., Knee Surg Sports Traumatol Arthrosc., 2012), of grade I and II according to their classification by Kellgren and Lawrence, in the absence of osteophytes.

Maioregen Chondro+ is indicated for the treatment of single and multiple chondral lesions, with no or slight alteration of the subchondral osseous tissue, of grade III and IV according to their classification by Outerbridge, whether traumatic, post-traumatic or degenerative.

Maioregen devices may be used only for cartilaginous or osteochondral defects of size less than 10 cm² (each lesion), in patients aged 15 to 60 years.

AVAILABLE FORMATS

Maioregen devices are available in three configurations, of various shapes and sizes:

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Pre-operative treatment

It is advisable to perform pre-operative X-ray and MRI examinations for appropriate evaluation of the lesion characteristics (location, size, subchondral bone and cartilage quality) and the joint involved.

The patient must undergo prophylactic antibiotic therapy prior to surgery. In patients allergic to specific antibiotics, a substitutive prophylaxis should be considered.

Intra-operative aspects

Read the user instructions carefully before proceeding with the surgery.

Check that the package is undamaged and unopened and that the product has not expired. The device is intended exclusively for surgical use under sterile conditions.

Apply a tourniquet to the proximal region of the patient's thigh. Expose the area of the lesion with an arthrotomic or mini-invasive arthrotomic incision. Open the package only in a sterile environment.

Remove the damaged tissue with an appropriate surgical device (e.g. surgical curette, osteotome or mini-invasive arthrotomy instruments).

To ensure the mechanical stability and effectiveness of the device, use an appropriate surgical instrument to prepare a seat with side walls perpendicular to the base of the lesion, with a flat floor.

In particular, the seat must be of the following depth from the intact articular surface:

- Around 8.0 mm for **Maioregen Prime**,
- Around 5.0 mm for **Maioregen Slim**,
- Around 2.5 - 3.0 mm for **Maioregen Chondro+**.

If you wish to proceed with the use of **Maioregen Chondro+** stimulate the marrow (micro-fractures, nano-fractures, perforations) before you position the device, to enable marrow blood to flow from the subchondral compartment.

If necessary, shape the device based on the geometry of the prepared housing using a surgical scalpel for the cartilaginous layer and surgical scissors for the remaining part. Handle the device with care while shaping it, using clean, sterile gloves.

Make sure that you shape the device so that it is not larger than the housing.

Place the device in the prepared housing making sure that its rough surface is in contact with the bone part and that the smooth cartilaginous surface faces upwards. Exercise light pressure by hand. Make sure that the edges of the device are in contact with the receiving tissues and that they perfectly adapt to the housing.

We recommend that you further stabilise the implant with fibrin glue applied:

- to the edges of the device prior to positioning and at the surface inter-

face between the device and the articular cartilage after positioning, in the case of **Maioregen Prime and Maioregen Slim**;

- at the interface between the device and the articular cartilage after positioning, for **Maioregen Chondro+**.

You can also stabilise **Maioregen Chondro+** with reabsorbable polymer pins, applied at the corners or close to the edge of the device, exclusively in the case of lesions of the trochlea or the front part of the condyles, as an alternative to the use of fibrin glue.

Release the tourniquet; the stabilisation of the device is facilitated by the bleeding of the subchondral bone. Check the stability of the device by flexing and extending the joint repeatedly.

Close the surgical site according to standard procedures.

Post-operative treatment

In accordance with surgical procedures, the patient must undergo a suitable postoperative antibiotic therapy. In patients allergic to specific antibiotics, a substitutive prophylaxis should be considered.

Inform the patient of the rehabilitation treatments and therapies necessary for a proper post-operative course. Do not immobilize the limb after surgery. It is advisable to periodically check the treated site by means of NMR.

CONTRAINDICATIONS

Do not use the device in case of:

- advanced osteoarthritis;
- systemic or localised infection;
- inflammatory pathologies, both systemic (e.g. rheumatoid arthritis) and local (e.g. sinovitis);
- disorders of the immune system;
- hypercalcemia;
- degenerative or vascular bone disease (e.g. osteonecrosis);
- coagulation disorders;
- metabolic disorders;
- ongoing antineoplastic chemotherapy;
- diabetes;
- obesity (Body Mass Index, BMI > 30);

- chondral or osteochondral lesions exceeding 10 cm²;
- coinciding lesions;
- age under 15 years and over 60 years;
- thyroid function abnormalities or complications;
- systemic conditions affecting wound healing;
- known allergy to equine collagen and calcium phosphate salts.

USE IN SPECIFIC SUBJECTS

The effects of MaioRegen devices are not known in case of:

- renal diseases;
- pregnancy;
- breast-feeding;
- ongoing radiotherapy;
- neoplasia;
- cardiovascular diseases.

WARNINGS

The use of MaioRegen medical devices is reserved solely for qualified medical staff specialised in orthopaedics.

It is inadvisable to apply the device by means of mosaicplasty technique. The factors listed below should be considered when selecting the device and/or the therapeutic approach:

- age of the patient;
- quality of the bone and cartilaginous tissues;
- localisation of the defect;
- loading conditions;
- presence of temporary local inflammation of the joint;
- joint alignment.

MaioRegen devices are single-use only; do not use any remains for other operations or other patients but dispose of it as contaminated hospital waste according to the regulations in force at the hospital facility.

The product is not to be used if its packaging is damaged. The devices cannot be resterilised. Do not use the product after its expiry date.

Open the package carefully to prevent damaging the device. Avoid the dispersion of device particles outside the implantation site. If necessary, shape the device to size, to avoid using multiple fragments in treating a single lesion.

The package contains some labels bearing the traceability data of the product; these should be applied on the patient's clinical record.

Should the device be removed during the intra-operative stage, do not reuse it. Use a new device. In the event of allergic reaction, remove the device and perform an appropriate therapeutic treatment.

Avoid arthroscopic procedures at the implant site during the month following the operation.

The treatment of the defect may lead to unsatisfactory results in terms of regeneration of new cartilaginous or osteochondral connective tissue: the outcome of the treatment depends on the individual response of the patient and may vary in relation to the patient's observance of the post-operational rehabilitative protocol.

PRECAUTIONS

When placing the device in the housing prepared in the tissue, make sure that its rough surface is in contact with the bone part and that the smooth cartilaginous surface faces upwards.

To avoid compromising the integrity of the product, shape it to size so that it is not compressed excessively once positioned on the defect.

Maioregen devices do not guarantee mechanical support in the first post-implant weeks; the physician must prescribe the use of appropriate mobility accessories.

SIDE EFFECTS

During the post-operative phase, some complications related to the type of intervention may occur, such as prolonged swelling of the treated area, hyperthermia, joint adhesions and/or fibrosis.

STERILISATION

All Maioregen devices are supplied sterile (sterilized by gamma radiation: 25kGy). Each device is a single-use product and cannot be re-sterilised. Its reuse, in whole or in part, may involve the risk of cross-contamination and the risk of infection at the implant site.

STORAGE

Store the product in a cool and dry place away from light and heat sources (+10 °C / + 30 °C).

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FIN-CERAMICA FAENZA S.p.A.