RegenOss Spine € 0373

Italian IT

English EN

French FR

German DE

Spanish ES

Portuguese PT

Turkish TR

Dutch NL

Bulgarian BG

Slovak SK

Czech CZ

Slovenian SLO

	<p>Fabbricante Manufacturer Fabricant Hersteller Fabricante Fabricante Üretici Fabrikant Производител Výrobca Výrobce Proizvajalec</p>
	<p>Non risterilizzare Do not resterilize Ne pas restériliser Nicht reesterilisieren No volver a esterilizar Não submeter à reesterilização Sterilize etmeyiniz Niet opnieuw steriliseren Да не се стерилизира повторно Nesterilizujte opakovane Ne re-sterilizovat Ne sterilizirajte znova</p>
	<p>Non riutilizzare Do not re-use Ne pas réutiliser Nicht wiederverwenden No volver a usar Não reutilizar Yeniden kullanmayınız Niet hergebruiken Да не се използва повторно Neroužívajte opakovane Pouze na jednorázové užití Samo za enkratno uporabo</p>
	<p>Non utilizzare se la confezione è danneggiata Do not use if the package is damaged Ne pas utiliser si l'emballage est endommagé Nicht benutzen, wenn die Verpackung beschädigt ist No se debe utilizar si el envase está dañado Não utilizar se a embalagem estiver danificada Ambalaj hasarlı ise kullanmayınız. Niet gebruiken als de verpakking beschadigd is Да не се употребява, ако опаковката е повредена! Neroužívajte, ak je obal poškodený Při poškozeném obalu nelze použít Izdelka ne uporabljajte, če je embalaža poškodovana</p>

	<p>Sterile, incluso il metodo di sterilizzazione (radiazioni ionizzanti) Sterile, including the method of sterilization (ionizing radiations) Stérile, incluant la méthode de stérilisation (radiations ionisantes) Steril, samt Sterilisationsmethode (ionisierende Strahlungen) Estéril, incluido el método de esterilización (radiaciones ionizantes) Estéril, incluindo o método de esterilização (radiações ionizantes) Steril, sterilizasyon yöntemi (iyonlaştırıcı radyasyon) dahil Steriel, inclusief sterilisatiemethode (ioniserende straling) Стерилен, включително метод на стерилизация (ионизиращо лъчение) Sterilné, metóda sterilizácie (ionizujúce žiarenie) Sterilní, včetně sterilizační metody (ionizující záření) Sterilno, vključno z načinom sterilizacije (ionizirajoče sevanje)</p>
	<p>Utilizzare entro Use by Date limite d'utilisation Mindestens haltbar bis Fecha de caducidad Usar antes de Son kullanma tarihi Te gebruiken tot Используй до Dátum expirácie Datum použitelnosti Rok uporabnosti</p>
	<p>Marcatura CE CE marking Marquage CE CE-Kennzeichnung Marcado CE Marcação CE CE işareti CE-markering CE маркировка CE značka notifikovaného organu CE značka Oznaka CE</p>
	<p>Numero di lotto Batch code Numéro de lot Chargennummer Número de lote Número de lote Parti numarası Lotnummer Партиден № Číslo šarže Číslo šarže Številka serije</p>

	<p>Numero di catalogo Catalogue number Numéro de catalogue Katalog-Nummer Número catálogo Número de catálogo Katalog numarası Catalogusnummer Каталоген № Katalógové číslo Katalógové číslo Kataloška številka</p>
	<p>Numero di pezzi Number of pieces Nombre de pièces Stückzahl Número de piezas Número de unidades Adet sayısı Aantal stuks Бройки от продукта Počet kusov v balení Počet kusů Število kosov</p>
	<p>Attenzione, vedere le istruzioni per l'uso Warning: see instructions for use Attention, voir les instructions pour l'emploi Achtung, siehe die Gebrauchshinweise Cuidado, consultar las instrucciones de uso Atenção, consultar as instruções de uso Dikkat, kullanım talimatlarına bakınız Let op, zie de gebruiksaanwijzing Внимание, прочетете инструкцията за употреба! Pozor, přeštudujte si priložený návod na použití Pozor, přečtěte si návod k použití Opozorilo: upoštevajte navodila za uporabo</p>
	<p>Consultare le istruzioni operative Consult the operating instructions Voir les instructions pour l'utilisation Bedienungsanleitung beachten Consultar las instrucciones operativas Consultar as instruções operacionais İşletim talimatlarına başvurunuz Lees de gebruiksaanwijzing Консултирайте се с инструкцията за употреба! Dodržujte návod na použití Prostudujte si pokyny Preverite navodila za uporabo</p>

RegenOss Spine

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DESCRIZIONE

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Bone substitutes

DESCRIPTION

RegenOss is a medical device for bone defect reconstruction.

RegenOss is an implantable, reabsorbable, non-active medical device, acting as bone substitute. It is composed of type-I collagen (of equine origin) and biomimetic hydroxyapatite whose chemical, chemical-physical and biological features are equivalent to those of human bone mineral constituent.

Collagen is a fibrous protein and represents the most important structural component of many human connective tissue extracellular matrixes, while hydroxyapatite is the main bone mineral element (65% by weight).

The presence of type-I collagen makes this device chemically and structurally biomimetic and ensures great in-site stability.

RegenOss is a flexible composite that can reproduce the bone anatomical structure. RegenOss, thanks to its special configuration and biomimetic chemical composition, can promote bone tissue regeneration processes.

INDICATIONS for USE

RegenOss is indicated for recovering bone substance loss, in particular:

- in presence of bone defects due to trauma, degenerative diseases of the skeletal system or surgical excisions, in combination with osteosynthesis implants and/or prosthesis;
- for filling defects due to autologous bone removal;
- during spinal arthrodesis procedures (interbody fusion and posterolateral fusion), and for spondylolisthesis;

The device acts as a bone substitute that is reabsorbed during the remodelling process of the newly-formed tissue.

RegenOss can be used in combination with saline fluids to improve its flexibility. The device can be used with surgical glues to increase stability in case of application on tissue defects with uneven or jagged surfaces.

AVAILABLE SIZES

The device is available in various pre-shaped forms and sizes.

INSTRUCTIONS for USE

Instructions for the surgical use of RegenOss are described below.

The following instructions are to be considered only as recommendations for the surgical use of RegenOss, they are not intended to modify or replace standard procedures for the treatment of bone defects. Surgical procedures involving the use of bone substitutes may have variable outcomes. The factors to be taken into account when choosing a bone substitute and the surgical technique are as follows: age and bone quality of the patient, location of the defect, expected loading condition, level and quality of tissue vascularity near the defect, possibility for the substitute to be in direct contact with the vital bone tissue adjacent to the defect, chances to ensure the complete filling of the bone defect, possibility to obtain correct and sufficient primary stabilisation of the bone substitute, chances to obtain complete coverage of the bone substitute.

Pre-operative treatment

It is advisable to carry out pre-operative radiographs to evaluate the conditions of the bone segment and the size of the defect to be treated as well as to choose the type and amount of product to apply and the fixation methods and devices. In the event the surgeon decides to proceed with surgery, the patient must undergo antibiotic prophylaxis. An alternative prophylaxis is essential to be used in patients allergic to penicillins.

It is important to carefully check the absence of infections or inflammations prior to surgery.

Intra-operative aspects

The device must be used only for surgeries performed in sterile operating environments, adopting every precaution. The package is to be opened after verifying its integrity and only in sterile environments.

RegenOss does not interfere with normal surgical procedures.

After careful cleaning of the bone surface of the defect to be treated, the product can be directly applied.

It is essential to prevent any dispersion of material outside the implant site but the amount of product used must be enough to fill the entire bone defect.

If necessary, shape the device according to the site by means of sterile forceps or scalpels. The quantity of device to be used must not exceed the size of the defect.

Place the device in the bone defect to be treated and stabilise it, pressing gently.

If necessary, use surgical glue for further stabilisation. Bleeding of the bone compartment may facilitate stabilisation of the device.

Excess material must be carefully removed making sure not to dislocate the bone substitute placed in the defect. Avoid irrigating the defect site once the substitute has been positioned.

Once the filling is completed, close the surgical site following the standard procedure.

RegenOss does not have sufficient mechanical resistance to bear direct loads in the immediate post-op phase. If applied in areas that require immediate loading, it is important to use suitable fixation means or synthetic systems.

Post-operative treatment

The patient should undergo a suitable post-operative antibiotic therapy according to the surgical procedures.

The surgeon must give the patient all the instructions related to the specific application to guarantee a proper post-operative course.

However, it is important to verify the newly formed bone tissue radiographically before subjecting the area to direct loading.

CONTRAINDICATIONS

Do not use the device in case of:

- systemic or localised infection;
- inflammatory or autoimmune disease;
- hypercalcemia;
- degenerative bone disease;
- coagulation disorders;
- metabolic disorders;
- diabetes;
- thyroid function abnormalities or complications;
- systemic conditions affecting wound healing;
- allergies to collagen of equine origin and to calcium phosphate salts.

The product must not be implanted in areas where complete coverage with soft tissue is not possible.

USE IN SPECIFIC SUBJECTS

RegenOss efficiency is unknown in case of:

- renal diseases;
- pregnancy;
- breast-feeding;
- ongoing radiotherapy;
- neoplasia;
- cardiovascular diseases.

No side effects due to the use RegenOss in presence of cancer have been

reported. However, patients undergoing bone filler treatment following tumour removal must be carefully monitored.

WARNINGS

RegenOss is intended to be used only by qualified and specialised medical personnel.

RegenOss does not have sufficient mechanical resistance to bear direct loads before regeneration of the newly formed tissue. If applied in areas that require immediate stabilisation and/or loading, it is important to use fixation means or synthetic systems.

Do not use the product after its expiry date.

RegenOss is a single-use product. Do not use any remains for other operations or other patients but dispose of it as contaminated hospital waste according to the regulations in force at the hospital facility. The product must not be re-sterilised.

The manufacturer is not liable for risks deriving from processes to recondition the device; reconditioning of the device implies the loss of its features and functions.

The product is not to be used if its packaging is damaged. The sterility of the product is guaranteed only if the package is unopened and undamaged.

The package contains an additional label displaying product tracking information; the implant must be tracked with the label, which must be applied to the patient's medical records.

In case of allergic reaction, the device must be removed and the patient must be treated with an adequate therapy.

The factors to be taken into account when choosing the device and the surgical technique are as follows:

- age of the patient;
- bone tissue quality;
- localisation of the defect;
- loading conditions.

Surgical procedures involving the use of bone substitutes may have variable outcomes, which also depend on the patient's compliance with the post-op protocol.

PRECAUTIONS

In order not to compromise the porosity of the product, do not use an excessive amount of material or press the device too much when placing it in the defect

site. RegenOss does not guarantee immediate post-implantation mechanical support, therefore it is important to use an adequate stabilisation system.

STERILIZATION

All RegenOss formats are supplied sterile (sterilized by gamma 25 kGy rays). The device is a single-use product and cannot be resterilised. Its reuse, in whole or in part, may involve the risk of cross-contamination and the risk of infection at the implant site.

STORAGE

Keep the product in a cool and dry place, protected from light and from heat sources (+10 °C / +30 °C).

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FIN-CERAMICA FAENZA S.p.A.