SINTlife Spine

CE 0373

Italian IT

English EN

French FR

German DE

Spanish ES

Portuguese PT

Turkish TR

Dutch NL

Bulgarian BG

Slovak SK

Czech CZ

Slovenian SLO

	<p>Fabbricante Manufacturer Fabricant Hersteller Fabricante Fabricante Üretici Fabrikant Производител Výrobca Výrobce Proizvajalec</p>
	<p>Non risterilizzare Do not resterilize Ne pas restériliser Nicht reesterilisieren No volver a esterilizar Não submeter à reesterilização Sterilize etmeyiniz Niet opnieuw steriliseren Да не се стерилизира повторно Nesterilizujte opakovane Ne re-sterilizovat Ne sterilizirajte znova</p>
	<p>Non riutilizzare Do not re-use Ne pas réutiliser Nicht wiederverwenden No volver a usar Não reutilizar Yeniden kullanmayınız Niet hergebruiken Да не се използва повторно Neroužívejte opakovane Pouze na jednorázové užití Samo za enkratno uporabo</p>
	<p>Non utilizzare se la confezione è danneggiata Do not use if the package is damaged Ne pas utiliser si l'emballage est endommagé Nicht benutzen, wenn die Verpackung beschädigt ist No se debe utilizar si el envase está dañado Não utilizar se a embalagem estiver danificada Ambalaj hasarlı ise kullanmayınız. Niet gebruiken als de verpakking beschadigd is Да не се употребява, ако опаковката е повредена! Neroužívejte, ak je obal poškozený Při poškozeném obalu nelze použít Izdelka ne uporabljajte, če je embalaža poškodovana</p>

	<p>Sterile, incluso il metodo di sterilizzazione (radiazioni ionizzanti) Sterile, including the method of sterilization (ionizing radiations) Stérile, incluant la méthode de stérilisation (radiations ionisantes) Steril, samt Sterilisationsmethode (ionisierende Strahlungen) Estéril, incluido el método de esterilización (radiaciones ionizantes) Estéril, incluindo o método de esterilização (radiações ionizantes) Steril, sterilizasyon yöntemleri (iyonlaştırıcı radyasyon) dahil Steriel, inclusief sterilisatiemethode (ioniserende straling) Стерилен, включително метод на стерилизация (ионизиращо лъчение) Sterilné, metóda sterilizácie (ionizujúce žiarenie) Sterilní, včetně sterilizační metody (ionizující záření) Sterilno, vključno z načinom sterilizacije (ionizirajoče sevanje)</p>
	<p>Utilizzare entro Use by Date limite d'utilisation Mindestens haltbar bis Fecha de caducidad Usar antes de Son kullanma tarihi Te gebruiken tot Используй до Dátum expirácie Datum použitelnosti Rok uporabnosti</p>
	<p>Marcatura CE CE marking Marquage CE CE-Kennzeichnung Marcado CE Marcação CE CE işareti CE-markering CE маркировка CE značka notifikovaného organu CE značka Oznaka CE</p>
	<p>Numero di lotto Batch code Numéro de lot Chargennummer Número de lote Número de lote Parti numarası Lotnummer Партиден № Číslo šarže Číslo šarže Številka serije</p>

	<p>Numero di catalogo Catalogue number Numéro de catalogue Katalog-Nummer Número catálogo Número de catálogo Katalog numarası Catalogusnummer Каталоген № Katalógové číslo Katalógové číslo Kataloška številka</p>
	<p>Numero di pezzi Number of pieces Nombre de pièces Stückzahl Número de piezas Número de unidades Adet sayısı Aantal stuks Бройки от продукта Počet kusov v balení Počet kusů Število kosov</p>
	<p>Attenzione, vedere le istruzioni per l'uso Warning: see instructions for use Attention, voir les instructions pour l'emploi Achtung, siehe die Gebrauchshinweise Cuidado, consultar las instrucciones de uso Atenção, consultar as instruções de uso Dikkat, kullanım talimatlarına bakınız Let op, zie de gebruiksaanwijzing Внимание, прочетете инструкцията за употреба! Pozor, přeštudujte si příložený návod na použití Pozor, přečtěte si návod k použití Opozorilo: upoštevajte navodila za uporabo</p>
	<p>Consultare le istruzioni operative Consult the operating instructions Voir les instructions pour l'utilisation Bedienungsanleitung beachten Consultar las instrucciones operativas Consultar as instruções operacionais İşletim talimatlarına başvurunuz Lees de gebruiksaanwijzing Консултирайте се с инструкцията за употреба! Dodržujte návod na použití Prostudujte si pokyny Preverite navodila za uporabo</p>

SINTlife Spine

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SINTlife Spine

Bone substitutes made of hydroxyapatite enriched with magnesium

DESCRIPTION

The mineral constituent of bone structure is composed of hydroxyapatite. Hydroxyapatite is found in human bones in percentages ranging between 65% and 80%.

SINTlife is a range of bone substitutes made of a new generation biomimetic hydroxyapatite enriched with magnesium ions in a similar amount as in human bones. Therefore, SINTlife chemical and chemical-physical features (crystallinity) are equivalent to those of human bone mineral constituent.

Thanks to their specific biomimetic chemical composition and their microstructural features, SINTlife medical devices promote the new bone formation process, demonstrating high osteointegration kinetics.

INDICATIONS for USE

SINTlife may be used to recover bone substance loss, in particular:

- in presence of bone defects due to trauma, degenerative diseases of the skeletal system or surgical excisions, in combination with osteosynthesis implants and/or prosthesis;
- for filling defects due to autologous bone removal or to tumour and cysts removal;
- during spinal arthrodesis procedures (interbody fusion and posterolateral fusion), and for spondylolisthesis;

SINTlife can also be used:

- in combination with other Finceramica bone substitutes;
- in combination with autologous, homologous and synthetic bone;
- mixed with peripheral blood, marrow blood, platelet gel (PRP) and growth factors;
- mixed with haemostatic agents and biological glues;
- mixed with antibiotics.

AVAILABLE SIZES

SINTlife is available in form of microgranules, injectable paste and malleable putty.

INSTRUCTIONS for USE

The following instructions are to be considered as recommendations for the surgical use of SINTLife, they are not intended to modify or replace standard procedures for the treatment of bone defects. Surgical procedures involving the use of bone substitutes may have variable outcomes. The factors to be taken into account when choosing a bone substitute and the surgical technique are as follows: age and bone quality of the patient, location of the defect, expected loading condition, level and quality of tissue vascularity near the defect, possibility for the substitute to be in direct contact with the living bone tissue adjacent to the defect, chances to ensure the complete filling of the bone defect, possibility to obtain correct and sufficient primary stabilisation of the bone substitute, chances to obtain complete coverage of the bone substitute.

Pre-operative treatment

It is advisable to carry out pre-operative radiographs to evaluate the conditions of the bone segment and the size of the defect to be treated as well as to choose the type and amount of product to apply and the fixation methods and devices. In the event the surgeon decides to proceed with surgery, the patient must undergo antibiotic prophylaxis. An alternative prophylaxis is essential to be used in patients allergic to penicillins.

It is important to carefully verify the absence of infections or inflammations prior to surgery.

Intra-operative aspects

Use SINTLife only for surgeries performed in sterile environments, adopting every precaution. The product must be opened only in sterile environments. SINTLife does not interfere with normal surgical procedures. After careful cleaning of the bone surface, the product can be directly applied onto the bone area to be treated.

It is essential to prevent any dispersion of material outside the implant site but the amount of product used must be enough to fill the entire bone defect. Excess material must be carefully removed making sure not to dislocate the bone substitute placed in the defect. Avoid irrigating the defect site once SINTLife has been placed.

Once the filling is completed, close the surgical site following the standard procedure.

SINTLife does not have sufficient mechanical resistance to bear direct loads in the immediate post-op phase. If applied in areas that require immediate loading, it is important to use suitable fixation means or synthetic systems.

Post-operative treatment

The patient should undergo a suitable post-operative antibiotic therapy according to the surgical procedures. The surgeon must give the patient all the instructions related to the specific application to guarantee a proper post-operative course.

However, it is important to verify the newly formed bone tissue radiographically before subjecting the area to direct loading.

CONTRAINDICATIONS

Do not use SINTlife in presence of inflammations and infections in the surgical site to be treated.

As for neurosurgical and spine applications, SINTlife must not be used in presence of even suspected dura mater injuries.

Moreover, it should not be used in patients that are hypersensitive to calcium phosphate salts or affected by: severe diabetes and demineralisation, serious impairment or liver and renal diseases, hypercalcaemia, thyroid and parathyroid complications, altered coagulation, infections due to previous surgeries, metabolic disorders of any kind, nor in patients undergoing steroid and anti-coagulant therapies.

The product must not be implanted in areas where complete coverage with soft tissue is not possible.

No side effects due to the use SINTlife in the presence of cancer have been reported. However, patients undergoing bone filling treatment following tumour removal should be carefully monitored.

WARNINGS

SINTlife does not have sufficient mechanical resistance to bear direct loads before regeneration of the newly formed tissue. If applied in areas that require immediate stabilisation and/or loading, it is important to use fixation means or synthetic systems.

SINTlife is intended to be used only by qualified and specialised medical personnel.

Do not use the product if expired.

The product is not to be used if its packaging is damaged.

SINTlife is a single-use product. Do not use any remains for other operations or other patients but dispose of it as contaminated hospital waste according to the regulations in force at the hospital facility. The product must not be re-sterilised. The manufacturer is not liable for risks deriving from processes to recondition the device; reconditioning of the device implies the loss of its features and functions.

The package contains an additional label displaying product tracking information; the implant must be tracked with the label, which must be applied to the patient's medical records.

SIDE EFFECTS

No side effects have been reported following the use of this product according to the correct indications.

STERILIZATION

All SINTlife formats are supplied sterile (sterilised by gamma rays, at a dose of 25 kGy). These are single-use products and must not be re-sterilised. Its reuse, in whole or in part, involves the risk of cross contamination and the danger of infection at the implant site.

STORAGE

Keep the product in a cool and dry place, protected from light and from heat sources (+10 °C / +40 °C).

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FIN-CERAMICA FAENZA S.p.A.